

THE PHOTOCONDUCTIVE STIMULATION DEVICE (PSD)

“Non-invasive, Controlled, Depolarization of Excitable Cells”



Advancing life science research through enabling rapid, non-invasive, targeted excitation of excitable cells quickly and easily

- Non-invasive photoconductive stimulation reduces excitotoxicity supporting extended periods of analysis
- Optical targeting allows excitation targeting from single cells to whole cell cultures
- Allows for the controlled variation of frequency and duration of stimulation
- The depolarization of a user-defined area of a dissociated culture, not just those cells positioned on or near an electrode
- Enables significant increase in experimental protocol flexibility and design by readily having more control over which cells are depolarized and when while maintaining cellular viability
- Depolarize cells with specific patterns of excitation, which allows for detailed investigation of excitable cells over time
- Significantly reduces the time and skill level required to obtain excitation, in that electrophysiology techniques and equipment are not required

Photoconductive Stimulation at a Glance:

- The photoconductive effect alters the conductivity of silicon with exposure to light
- Illuminate as small an area as approximately ~50 microns
- Stimulation protocols can be from 0.1 Hz up to 100 Hz and approximately up to a 17 volt pulse across entire surface of wafer
- Only illuminated neuron(s) depolarize on-demand
- Study cellular response using standard biochemical protocols or live imaging
- Non-invasive, preserves cytoskeleton, physiologically benign
- Concurrent electrophysiological recordings demonstrate high frequency stimulation possible
- Stimulation can result in controlled synaptic and network remodeling
- Networks formed have complex activity patterns
- Ion channel analysis
- Long term potentiation in the presence of specific controlled brain activity
- Depolarize cells in the presence of pharmaceutical compounds to screen for efficacy
- Long term analysis of cultures undergoing controlled depolarization
- Appropriate for use with any excitable cells, including those from the hippocampus, cortex and striatum

More Information:

Website: www.neurosilicon.com
E-mail: contact@neurosilicon.com



Neurosilicon

Question: Does the Photoconductive Stimulation Device (PSD) sense cellular activity?

Answer: No, the PSD does not sense any cellular activity. The device is used solely for depolarizing cells in a controlled manner. In some cases, sensing may be necessary. The PSD does operate alongside standard fluorescent imaging systems for this purpose.

Question: Do I need any “special” cells to use the PSD?

Answer: No, the PSD operates with any standard mammalian excitable cells, including neurons and cardiac myocytes.

Question: What kind of equipment do I need to prepare cells for use with the PSD?

Answer: The standard tools to culture cells are necessary. No special tools are required for this purpose. Neurosilicon does provide the silicon materials on which to culture the cells.

Question: How does the PSD stimulate cells without damaging or physically touching them?

Answer: The PSD uses a phenomenon called the Photoconductive Stimulation Effect. It is based in the principle that when visible light strikes silicon in combination with a specific current, a local field potential is generated. This LFP is what causes the cell to depolarize noninvasively without causing any damage to the cell.

Question: Do I need a microscope to use the PSD with?

Answer: Yes, the PSD does require a microscope for operation. The PSD only supports the use of a standard upright microscope with the light source coming from above.

Question: Do I need any special training to use the PSD?

Answer: No, the PSD is a standard neuroscience research tool that any skilled technician can use with the use of the appropriate documentation.

In Press:

- *Activation of microglia by neuronal activity: Results from a new in vitro paradigm based on neuronal-silicon interfacing technology.* **Brain Behav Immun**, July. (2009)
- *Astrocytic Ca²⁺ Waves Guide CNS Growth Cones to Remote Regions of Neuronal Activity.* **PLoS One**; 3(11): e3692. (2008)
- *Photoconductive Stimulation Of Neurons Cultured On Silicon Wafers.* **Nature Protocols** 1, 461 – 467. (2006)
- *Neuronal Networks And Synaptic Plasticity: Understanding Complex System Dynamics By Interfacing Neurons With Silicon Technologies.* **Journal of Experimental Biology**, 209 (Pt 12): 2312-9. (2006)
- *Remodeling Of Synaptic Actin Induced By Photoconductive Stimulation.* **Cell**, 107(5): 605-616. (2001)

More Information:

Website: www.neurosilicon.com
E-mail: contact@neurosilicon.com